

Nominal Power	450W	900W
kV Range	40~90kV	40~90kV
mA Range	2~5mA	2~10mA
Focal Spot Size	0.8mm	0.8mm
Anode Disk Target Angle	19°	19°
Communication	RS232/WIFI	RS232/WIFI
Operation Temperature	-20℃ ~ 40℃	-20℃ ~ 40℃
Dimensions (L*W*H)	266.4*218.4*157.7mm	266.4*218.4*157.7mm
Weight	≤2.8kg	≤3kg